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#### HPLC conditions

The samples were injected (5  $\mu$ L) into UHPLC connected to a photodiode array detector (Dionex Ultimate 3000), with a reverse-phase column (Ace Excel, 100\*3 mm, 1.7  $\mu$ m). The mobile phase consisted of (A) DDW with 1mM ammonium acetate and (B) Methanol containing 1mM ammonium acetate. The gradient was started with 40% B and kept isocratic for 2 min, then increased to 80% B in 8 min, then increased to 90% B in 2 minutes, and then kept isocratic for 5 min before returning to 40% B in 1 min and column was allowed to equilibrate at 40% B for 4 min before the next injection. The flow rate was 0.4 mL/min. column temperature was 35C.

#### LC/MS analysis

MS<sup>1</sup> and MS<sup>2</sup> analysis performed with Heated Electrospray ionization (HESI-II) source connected to a Q Exactive™ Plus Hybrid Quadrupole-Orbitrap™ Mass Spectrometer Thermo Scientific™. ESI capillary voltage was set to 3500 V, capillary temperature to 350C°, gas temperature to 250C° and Nitrogen gas flow to 35 mL/min. The mass spectra (m/z 150–900) were acquired using negative ion mode.

#### Data preprocessing

Peak determination, peak area integration, removing blank peaks, compound identification applied using Qn™ software.

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